







Fisher & Paykel Brevida Nasal Pillows Mask

Cleaning/Replacement Guidelines

	ResMed Product#	HCPCS Code	Replacement Schedule
 Brevida Nasal Pillows Mask (with headgear)	Brevida mask with headgear with: <ul style="list-style-type: none"> • XS-S pillow - BRE1SA • M-L pillow - BRE1MA • Fit Pack, XS-S & M-L pillows - BRE1SMA 	A7034 + A7035 (Mask + Headgear)	
 Brevida Nasal Pillows Mask	Brevida mask with: XS/S pillow - 400BRE131 M/L pillow - 400BRE132	A7034 (Mask)	Replace every <u>3 months</u> .
 Brevida Headgear	Brevida Headgear <ul style="list-style-type: none"> • 400BRE120 	A7035 (Headgear)	Replace every <u>6 months</u> .
 Brevida AirPillow	Brevida AirPillow size: <ul style="list-style-type: none"> • XS-S - 400BRE113 • M-L - 400BRE114 	A7033 (Pillows)	Replace at least every <u>3 months</u> . May replace every 2 weeks if needed.
Cleaning Schedule	Wash <u>weekly</u> with warm water and mild soap and air dry. May rinse pillow, mask, elbow, swivel, and tube daily if needed.		

DME Provider:
Phone Number: