







Philips DreamWear Silicone Pillows Mask

Cleaning/Replacement Guidelines

	ResMed Product#	HCPCS Code	Replacement Schedule
 DreamWear Silicone Pillows	<p>SMALL frame with: S cushion - 1146410, M cushion - 1146411 L cushion - 1146412, MW cushion - 1146413</p> <p>MEDIUM frame with: S cushion - 1146414, M cushion - 1146415 L cushion - 1146416, MW cushion - 1146417 Fitpack with S, M, L, & MW cushions - 1146468</p> <p>LARGE frame with: S cushion - 1146448, M cushion - 1146449 L cushion - 1146450, MW cushion - 1146451</p>	A7034 + A7035 (Mask + Headgear)	
 DreamWear Silicone Pillows	<p>SMALL frame with: S cushion - 1116710, M cushion - 1116711 L cushion - 1116712, MW cushion - 1116713</p> <p>MEDIUM frame with: S cushion - 1116705, M cushion - 1116706 L cushion - 1116707, MW cushion - 1116708 Fitpack with S, M, L, & MW cushions - 1116701</p> <p>LARGE frame with: S cushion - 1116715, M cushion - 1116716 L cushion - 1116717, MW cushion - 1116718</p>	A7034 (Mask)	Replace every 3 months.
 DreamWear Headgear or DreamWear Headgear with Arms	<p>DreamWear Headgear</p> <ul style="list-style-type: none"> • 1116750 <p>DreamWear Headgear with Arms</p> <ul style="list-style-type: none"> • 1116751 	A7035 (Headgear)	Replace every 6 months.
 DreamWear Silicone Pillows Cushion	<p>Nuance Gel Pillows Cushion size:</p> <ul style="list-style-type: none"> • S - 1146472 • M - 1146473 • L - 1146474 • MW - 1146475 	A7033 (Pillows)	Replace at least every 3 months. May replace every 2 weeks if needed.
Cleaning Schedule	Wash weekly with warm water and mild soap and air dry. May rinse pillows, frame, elbow, and swivel daily if needed.		

DME Provider:
Phone Number: