







ResMed AirFit N20/N20 for Her Nasal Mask Cleaning/Replacement Guidelines

	ResMed Product#	HCPCS Code	Replacement Schedule
 AirFit N20 or AirFit N20 for Her Mask (with headgear)	AirFit N20 mask, standard headgear with cushion size: <ul style="list-style-type: none"> • S - 63503 • M - 63501 • L - 63502 AirFit N20 for Her mask, small headgear with cushion size: <ul style="list-style-type: none"> • S - 63500 	A7034 + A7035 (Mask + Headgear)	
 AirFit N20 Frame System	AirFit N20 frame with cushion size: <ul style="list-style-type: none"> • S - 63553 • M - 63555 • L - 63556 AirFit N20 for Her frame with cushion size: <ul style="list-style-type: none"> • S - 63554 	A7034 (Mask)	Replace every <u>3 months.</u>
 AirFit N20 Headgear	AirFit N20 Headgear: <ul style="list-style-type: none"> • Small - 63560 • Standard - 63561 • Large - 63562 AirFit N20 for Her Headgear: <ul style="list-style-type: none"> • Small - 63558 	A7035 (Headgear)	Replace every <u>6 months.</u>
 AirFit N20 Cushion	AirFit N20 cushion size: <ul style="list-style-type: none"> • S - 63550 • M - 63551 • L - 63552 	A7033 (Cushions)	Replace at least every <u>3 months.</u> May replace every 2 weeks if needed.
Cleaning Schedule	Wash <u>weekly</u> with warm water and mild soap and air dry. May rinse cushion, vent, and short tube daily if needed.		

DME Provider:
Phone Number: