







ResMed AirFit P30i Pillows Mask

Cleaning/Replacement Guidelines

	ResMed Product#	HCPCS Code	Replacement Schedule
 AirFit P30i Standard Starter Pack or AirFit P30i Small Starter Pack	AirFit P30i, standard frame (S, M, L pillows included) <ul style="list-style-type: none"> • 63850 AirFit P30i, small frame (S, M pillows included) <ul style="list-style-type: none"> • 63851 	A7034 + A7035 (Mask + Headgear)	
 AirFit P30i Frame System	Standard frame with: <ul style="list-style-type: none"> • S pillow - 63856 • M pillow - 63584 • L pillow - 63852 Small frame with: <ul style="list-style-type: none"> • S pillow - 63857 • M pillow - 63855 • L pillow - 63853 	A7034 (Mask)	Replace every <u>3 months.</u>
 AirFit P30i Headgear	AirFit P30i Headgear <ul style="list-style-type: none"> • 63814 	A7035 (Headgear)	Replace every <u>6 months.</u>
 AirFit P30i Pillow Cushion	AirFit P30i pillow size: <ul style="list-style-type: none"> • S - 62931 • M - 62932 • L - 62933 	A7033 (Pillows)	Replace at least every <u>3 months.</u> May replace every 2 weeks if needed.
Cleaning Schedule	Wash <u>weekly</u> with warm water and mild soap and air dry. May rinse pillows daily if needed.		

DME Provider:
Phone Number: