







## Philips DreamWisp Nasal Mask

### Cleaning/Replacement Guidelines

	Philips Product#	HCPCS Code	Replacement Schedule
 <p><b>DreamWisp mask (with headgear)</b></p>	Medium connector with cushion size: <ul style="list-style-type: none"> <li>• P - 1137932</li> <li>• S - 1137933</li> <li>• M - 1137934</li> <li>• L - 1137935</li> <li>• XL - 1137936</li> <li>• Fitpack with S, M, L - 1137916</li> </ul>	<b>A7034 + A7035</b> (Mask + Headgear)	
 <p><b>DreamWisp mask</b></p>	Medium connector with cushion size: <ul style="list-style-type: none"> <li>• P - 1137922</li> <li>• S - 1137923</li> <li>• M - 1137924</li> <li>• L - 1137925</li> <li>• XL - 1137926</li> <li>• Fitpack with S, M, L - 1137918</li> </ul>	<b>A7034</b> (Mask)	Replace every <b>3 months.</b>
 <p><b>DreamWisp Headgear</b></p>	DreamWisp Headgear: <ul style="list-style-type: none"> <li>• Reduced size - 1137966</li> <li>• Standard - 1137964</li> <li>• Large - 1137965</li> </ul>	<b>A7035</b> (Headgear)	Replace every <b>6 months.</b>
 <p><b>DreamWisp Cushion</b></p>	DreamWisp cushion size: <ul style="list-style-type: none"> <li>• P - 1137956</li> <li>• S - 1137957</li> <li>• M - 1137958</li> <li>• L - 1137959</li> <li>• XL - 1137960</li> </ul>	<b>A7033</b> (Cushions)	Replace at least every <b>3 months.</b>  May replace every 2 weeks if needed.
<b>Cleaning Schedule</b>	Wash <b>weekly</b> with warm water and mild soap and air dry. May wash non-fabric parts daily if needed.		

DME Provider:  
Phone Number: