



## ResMed AirTouch F20 Full Face Mask

### Cleaning/Replacement Guidelines

	ResMed Product#	HCPCS Code	Replacement Schedule
 <b>AirTouch F20</b> or <b>AirTouch F20 for Her mask</b> (with headgear)	AirTouch F20 mask, <b>standard</b> headgear with cushion size: <ul style="list-style-type: none"> <li>• S - <b>63000</b></li> <li>• M - <b>63001</b></li> <li>• L - <b>63002</b></li> </ul> AirTouch F20 for Her mask, <b>small</b> headgear with cushion size: <ul style="list-style-type: none"> <li>• S - <b>63003</b></li> <li>• M - <b>63004</b></li> </ul>	<b>A7034 + A7035</b> (Mask + Headgear)	
 <b>AirTouch F20 Frame System</b>	AirTouch F20 frame with cushion size: <ul style="list-style-type: none"> <li>• S - <b>63021</b></li> <li>• M - <b>63022</b></li> <li>• L - <b>63023</b></li> </ul> AirTouch F20 for Her frame with cushion size: <ul style="list-style-type: none"> <li>• S - <b>63024</b></li> <li>• M - <b>63025</b></li> </ul>	<b>A7034</b> (Mask)	Replace every <b>3 months</b> .
 <b>AirTouch F20 Headgear</b>	AirTouch F20 Headgear <ul style="list-style-type: none"> <li>• Small - <b>63470</b></li> <li>• Standard - <b>63471</b></li> <li>• Large - <b>63472</b></li> </ul> AirTouch F20 for Her Headgear <ul style="list-style-type: none"> <li>• Small - <b>63473</b></li> </ul>	<b>A7035</b> (Headgear)	Replace every <b>6 months</b> .
 <b>AirTouch F20 Cushion</b>	AirTouch F20 cushion size: <ul style="list-style-type: none"> <li>• S - <b>63028</b></li> <li>• M - <b>63029</b></li> <li>• L - <b>63030</b></li> </ul>	<b>A7033</b> (Cushion)	Replace every <b>month</b> .
<b>Cleaning Schedule</b>	<b>DO NOT WASH AIRTOUCH F20 CUSHION.</b> Wash frame and headgear <u>weekly</u> with warm water and mild soap and air dry. May wipe cushion with moist wipes daily.		

**DME Provider:**  
**Phone Number:**