







Philips DreamWear Full Face Mask

Cleaning/Replacement Guidelines

	Philips Product#	HCPCS Code	Replacement Schedule
 DreamWear Full Face Mask (with headgear)	DreamWear Full Face Mask Fit Pack, MEDIUM frame with S, M, L, MW cushions - 1133400 DreamWear Full Face Mask, SMALL & MEDIUM frame with cushion size: S - 1133375 , M - 1133376 , L - 1133377 , MW - 1133378 DreamWear Full Face Mask, SMALL frame with cushion size: S - 1133385 , M - 1133386 , L - 1133387 , MW - 1133388 DreamWear Full Face Mask, MEDIUM frame with cushion size: S - 1133380 , M – 1133381 , L - 1133382 , MW - 1133383 DreamWear Full Face Mask, LARGE frame with cushion size: S - 1133390 , M – 1133391 , L - 1133392 , MW - 1133393	A7034 + A7035 (Mask + Headgear)	
 DreamWear Full Face Mask	DreamWear mask, SMALL frame with cushion size: S - 1133410 , M – 1133411 , L - 1133412 MW - 1133413 DreamWear mask, MEDIUM frame with cushion size: S - 1133405 , M – 1133406 , L - 1133407 , MW - 1133408 DreamWear mask, LARGE frame with cushion size: S - 1133415 , M – 1133416 , L - 1133417 , MW - 1133418	A7034 (Mask)	Replace every 3 months .
 DreamWear Full Face Headgear	DreamWear Full Face Headgear • 1133450	A7035 (Headgear)	Replace every 6 months .
 DreamWear Full Face Cushion	DreamWear Full Face Cushion size: • S - 1133430 • M - 1133431 • L - 1133432 • MW - 1133433	A7033 (Pillows)	Replace at least every 3 months . May replace every month if needed.
Cleaning Schedule	Wash weekly with warm water and mild soap and air dry. May wash non-fabric parts daily if needed.		

DME Provider:
Phone Number: