







ResMed AirFit N30 Nasal Mask

Cleaning/Replacement Guidelines

	ResMed Product#	HCPCS Code	Replacement Schedule
 AirFit N30 (with headgear)	AirFit N30 mask with cushion size: <ul style="list-style-type: none"> • S - 64222 • M - 64223 • SW - 64224 AirFit N30 Starter pack with S, M, SW: <ul style="list-style-type: none"> • 64201 	A7034 + A7035 (Mask + Headgear)	
 AirFit N30 Frame System	AirFit N30 frame with cushion size: <ul style="list-style-type: none"> • S - 64210 • M - 64211 • SW - 64212 	A7034 (Mask)	Replace every 3 months .
 AirFit N30 Headgear	AirFit N30 Headgear: <ul style="list-style-type: none"> • 64216 	A7035 (Headgear)	Replace every 6 months .
 AirFit N30 Cushion	AirFit N30 cushion size: <ul style="list-style-type: none"> • S - 64213 • M - 64214 • SW - 64215 	A7033 (Cushion)	Replace at least every 3 months . May replace every 2 weeks if needed.
Cleaning Schedule	Wash weekly with warm water and mild soap and air dry. May rinse cushions daily if needed.		

DME Provider:
Phone Number: