







ResMed AirFit N30i Nasal Mask

Cleaning/Replacement Guidelines

	ResMed Product#	HCPCS Code	Replacement Schedule
 <p>AirFit N30i mask (with headgear)</p>	<p>AirFit N30i Standard Starter pack (Standard frame with S, M, SW cushions)</p> <ul style="list-style-type: none"> 63800 <p>AirFit N30i, Small Starter pack (Small frame with S, M, SW cushions)</p> <ul style="list-style-type: none"> 63801 	<p>A7034 + A7035 (Mask + Headgear)</p>	
 <p>AirFit N30i Frame System</p>	<p>AirFit N30i, standard frame with cushion size:</p> <ul style="list-style-type: none"> S - 63802 M - 63804 W - 63806 SW - 63808 <p>AirFit N30i, Small frame with cushion size:</p> <ul style="list-style-type: none"> S - 63803 M - 63805 W - 63807 SW - 63809 	<p>A7034 (Mask)</p>	<p>Replace every 3 months.</p>
 <p>AirFit N30i Headgear</p>	<p>AirFit N30i Headgear</p> <ul style="list-style-type: none"> 63814 	<p>A7035 (Headgear)</p>	<p>Replace every 6 months.</p>
 <p>AirFit N30i Cushion</p>	<p>AirFit N30i cushion size:</p> <ul style="list-style-type: none"> S - 63813 M - 63811 W - 63812 SW - 63810 	<p>A7033 (Cushion)</p>	<p>Replace at least every 3 months.</p> <p>May replace every 2 weeks if needed.</p>
<p>Cleaning Schedule</p>	<p>Wash weekly with warm water and mild soap and air dry. May rinse cushions daily if needed.</p>		

DME Provider:
Phone Number: