







Philips Wisp Nasal Mask

Cleaning/Replacement Guidelines

	Philips Product#	HCPCS Code	Replacement Schedule
 Wisp mask (with headgear)	<p>Wisp CLEAR frame with cushion size: P - 1118060, S/M - 1118061 L - 1118062, XL - 1118063 Fitpack with S/M, L, XL cushions - 1094050</p> <p>Wisp FABRIC frame with cushion size: P - 1118064, S/M - 1118065 L - 1118066, XL - 1118067 Fitpack with S/M, L, XL cushions - 1094051</p>	A7034 + A7035 (Mask + Headgear)	
 Wisp mask	<p>Wisp CLEAR frame with cushion size: P - 1101550, S/M - 1101551 L - 1101552, XL - 1113957 Fitpack with S/M, L, XL cushions - 1094055</p> <p>Wisp FABRIC frame with cushion size: P - 1101553, S/M - 1101554 L - 1101555, XL - 1113958 Fitpack with S/M, L, XL cushions - 1094056</p>	A7034 (Mask)	Replace every <u>3 months.</u>
 Wisp Headgear	<p>Wisp Headgear:</p> <ul style="list-style-type: none"> • Reduced size - 1109307 • Standard - 1094082 • Large - 1094078 	A7035 (Headgear)	Replace every <u>6 months.</u>
 Wisp Cushion	<p>Wisp cushion size:</p> <ul style="list-style-type: none"> • P - 1094086 • S/M - 1094087 • L - 1094088 • XL - 1112031 	A7033 (Cushions)	Replace at least every <u>3 months.</u> May replace every 2 weeks if needed.
Cleaning Schedule	Wash <u>weekly</u> with warm water and mild soap and air dry. May wash non-fabric parts daily if needed.		

DME Provider:
Phone Number: