







ResMed AirFit F20/F20 for Her Full Face Mask

Cleaning/Replacement Guidelines

	ResMed Product#	HCPCS Code	Replacement Schedule
 <p>AirFit F20 or AirFit F20 for Her mask (with headgear)</p>	<p>AirFit F20 mask, standard headgear with cushion size:</p> <ul style="list-style-type: none"> • S - 63400 • M - 63401 • L - 63402 <p>AirFit F20 for Her mask, small headgear with cushion size:</p> <ul style="list-style-type: none"> • S - 63403 • M - 63404 	A7034 + A7035 (Mask + Headgear)	
 <p>AirFit F20 Frame System</p>	<p>AirFit F20 frame with cushion size:</p> <ul style="list-style-type: none"> • S - 63460 • M - 63461 • L - 63462 <p>AirFit F20 for Her frame with cushion size:</p> <ul style="list-style-type: none"> • S - 63463 • M - 63464 	A7034 (Mask)	Replace every 3 months .
 <p>AirFit F20 Headgear</p>	<p>AirFit F20 Headgear:</p> <ul style="list-style-type: none"> • Small - 63470 • Standard - 63471 • Large - 63472 <p>AirFit F20 for Her Headgear:</p> <ul style="list-style-type: none"> • Small - 63473 	A7035 (Headgear)	Replace every 6 months .
 <p>AirFit F20 Cushion</p>	<p>AirFit F20 cushion size:</p> <ul style="list-style-type: none"> • S - 63467 • M - 63468 • L - 63469 	A7033 (Cushion)	Replace at least every 3 months . May replace every month if needed.
Cleaning Schedule	Wash weekly with warm water and mild soap and air dry. May rinse cushion and frame daily if needed.		

DME Provider:
Phone Number: